


# Rho

BEST For ROSIN  
& LIVE RESIN

0.5ml 1.0ml




Pure Taste



**PURECORE®**  
Innovational  
Technology



**Highly**  
Customization



**POSTLESS**  
Cotton free  
Glue free



Window  
**Customized**  
0.5/1.0ml



**Anti-clog**  
Structure

## Features

- Special For Rosin and Live Resin;
- Pure And Mellow True To The Original Flavor;
- Highly Customizable, Create Your Exclusive Product;
- Anti-clog Structure helps to save your Extracts;
- Metal Housing;
- Inhale Activated;
- Experience The Excellence;
- “CLACK” To Save;

## Specification

| Product Name | Rho                   |
|--------------|-----------------------|
| Resistance   | 1.85Ω                 |
| Battery      | 265mAh                |
| Capacity     | 0.5ml/1.0ml           |
| Activation   | Inhale activation     |
| Charging     | Type-C                |
| Size         | 95.5x20x10mm          |
| Material     | PCTG, Stainless Steel |

## Color Plot



## Oil Filling Instruction

1. Insert your needle into your Rho oil tank through the top holes and hold upright while slowly filling. The maximum fill volume is 1.0ml/0.5ml .
  - ◆ Do not overfill or allow your oil to enter the airways on either side of the tanks.
2. After filling, immediately cap your Rho with the mouthpiece cover until it is tightly fitted.
  - ◆ Failing to attach the mouthpiece cover within 2 minutes may cause leakage.

